

FAM3C Polyclonal Antibody

Cat No: HR1AP12062

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | FAM3C Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 24kDa |
| GeneID?Human? | FAM3C ILE1 GS3786 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | family with sequence similarity 3 member C(FAM3C) Homo sapiens This gene is a member of the family with sequence similarity 3 (FAM3) family and encodes a secreted protein with a GG domain. A change in expression of this protein has been noted in pancreatic cancer-derived cells. [provided by RefSeq, Mar 2010], |