

CYR61 Polyclonal Antibody

Cat No: HR1AP12060

For research use only

Overview

Product Name	CYR61 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	CYR61 CCN1 GIG1 IGFBP10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cysteine rich angiogenic inducer 61(CYR61) Homo sapiens The secreted protein encoded by this gene is growth factor-inducible and promotes the adhesion of endothelial cells. The encoded protein interacts with several integrins and with heparan sulfate proteoglycan. This protein also plays a role in cell proliferation, differentiation, angiogenesis, apoptosis, and extracellular matrix formation. [provided by RefSeq, Sep 2011],</p>