

SFRP1 Polyclonal Antibody

Cat No: HR1AP12056

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | SFRP1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 34kDa |
| GeneID?Human? | SFRP1 FRP FRP1 SARP2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | secreted frizzled related protein 1(SFRP1) Homo sapiens This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. Members of this family act as soluble modulators of Wnt signaling; epigenetic silencing of SFRP genes leads to deregulated activation of the Wnt-pathway which is associated with cancer. This gene may also be involved in determining the polarity of photoreceptor cells in the retina. [provided by RefSeq, Sep 2009], |