

OLIG3 Polyclonal Antibody

Cat No: HR1AP12053

For research use only

Overview

Product Name	OLIG3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	OLIG3 BHLHB7 BHLHE20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,