



## SFRP4 Polyclonal Antibody

Cat No: HR1AP12052

For research use only

### Overview

Product Name	SFRP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	SFRP4 FRPHE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	secreted frizzled related protein 4(SFRP4) Homo sapiens Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression. [provided by RefSeq, Jul 2008],