

# ISLR2 Polyclonal Antibody

Cat No: HR1AP12048

For research use only

## Overview

Product Name	ISLR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	ISLR2 KIAA1465 LINX UNQ1885/PRO4329
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,similarity:Contains 5 LRR (leucine-rich) repeats.,