

PGCB Polyclonal Antibody

Cat No: HR1AP12047

For research use only

Overview

Product Name	PGCB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	BCAN BEHAB CSPG7 UNQ2525/PRO6018
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>brevican(BCAN) Homo sapiens This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011],</p>