

# PROM2 Polyclonal Antibody

Cat No: HR1AP12046

For research use only

## Overview

Product Name	PROM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	PROM2 PROML2 UNQ2521/PRO6014
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>prominin 2(PROM2) Homo sapiens This gene encodes a member of the prominin family of pentaspan membrane glycoproteins. The encoded protein localizes to basal epithelial cells and may be involved in the organization of plasma membrane microdomains. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009],</p>