

## TEFF2 Polyclonal Antibody

Cat No: HR1AP12045

For research use only

### Overview

Product Name	TEFF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	TMEFF2 HPP1 TENB2 TPEF UNQ178/PRO204
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>transmembrane protein with EGF like and two follistatin like domains 2(TMEFF2) Homo sapiens This gene encodes a member of the tomoregulin family of transmembrane proteins. This protein has been shown to function as both an oncogene and a tumor suppressor depending on the cellular context and may regulate prostate cancer cell invasion. Multiple soluble forms of this protein have been identified that arise from both an alternative splice variant and ectodomain shedding. Additionally, this gene has been found to be hypermethylated in multiple cancer types. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2015],</p>