

## NR0B1 Polyclonal Antibody

Cat No: HR1AP12035

For research use only

### Overview

Product Name	NR0B1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	NR0B1 AHC DAX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear receptor subfamily 0 group B member 1(NR0B1) Homo sapiens This gene encodes a protein that contains a DNA-binding domain. The encoded protein acts as a dominant-negative regulator of transcription which is mediated by the retinoic acid receptor. This protein also functions as an anti-testis gene by acting antagonistically to Sry. Mutations in this gene result in both X-linked congenital adrenal hypoplasia and hypogonadotropic hypogonadism. [provided by RefSeq, Jul 2008].