

RALA Polyclonal Antibody

Cat No: HR1AP12036

For research use only

Overview

Product Name	RALA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RALA RAL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RAS like proto-oncogene A(RALA) Homo sapiens The product of this gene belongs to the small GTPase superfamily, Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors. This gene encodes a low molecular mass ras-like GTP-binding protein that shares about 50% similarity with other ras proteins. [provided by RefSeq, Jul 2008],</p>