

# SALL1 Polyclonal Antibody

Cat No: HR1AP12034

For research use only

## Overview

Product Name	SALL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	145kDa
GeneID?Human?	SALL1 SAL1 ZNF794
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>spalt like transcription factor 1(SALL1) Homo sapiens The protein encoded by this gene is a zinc finger transcriptional repressor and may be part of the NuRD histone deacetylase complex (HDAC). Defects in this gene are a cause of Townes-Brocks syndrome (TBS) as well as bronchio-oto-renal syndrome (BOR). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>