

# DLL1 Polyclonal Antibody

Cat No: HR1AP12032

For research use only

## Overview

Product Name	DLL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	DLL1 UNQ146/PRO172
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	delta like canonical Notch ligand 1(DLL1) Homo sapiens DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis. It may play a role in cell-to-cell communication. [provided by RefSeq, Jul 2008],