

LRIG3 Polyclonal Antibody

Cat No: HR1AP12028

For research use only

Overview

Product Name	LRIG3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	123kDa
GeneID?Human?	LRIG3 LIG3 UNQ287/PRO326/PRO335
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 17 LRR (leucine-rich) repeats.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Widely expressed.,