

SPIT1 Polyclonal Antibody

Cat No: HR1AP12025

For research use only

Overview

Product Name	SPIT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	SPINT1 HAI1 UNQ223/PRO256
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serine peptidase inhibitor, Kunitz type 1(SPINT1) Homo sapiens The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008],