

FAM3B Polyclonal Antibody

Cat No: HR1AP12019

For research use only

Overview

Product Name	FAM3B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	FAM3B C21orf11 C21orf76 PRED44 UNQ320/PRO365
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Has delayed effects on beta-cell function, inhibiting basal insulin secretion from a beta-cell line in a dose-dependent manner.,miscellaneous:2 N-termini have been observed in the mature protein: the first at Glu-30, resulting from signal peptide cleavage, the second at Ser-46.,similarity:Belongs to the FAM3 family.,subcellular location:Localized in discrete vesicular and perinuclear structure.,tissue specificity:Highly expressed in the pancreas. Also found in the colon, kidney, prostate, small intestine and testis.,</p>