

XPP1 Polyclonal Antibody

Cat No: HR1AP12017

For research use only

Overview

Product Name	XPP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	XPNPEP1 XPNPEPL XPNPEPL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	X-prolyl aminopeptidase 1(XPNPEP1) Homo sapiens This gene encodes the cytosolic form of a metalloaminopeptidase that catalyzes the cleavage of the N-terminal amino acid adjacent to a proline residue. The gene product may play a role in degradation and maturation of tachykinins, neuropeptides, and peptide hormones. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Nov 2009],