

## CR3L1 Polyclonal Antibody

Cat No: HR1AP12015

For research use only

### Overview

Product Name	CR3L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	CREB3L1 OASIS PSEC0238
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cAMP responsive element binding protein 3 like 1(CREB3L1) Homo sapiens The protein encoded by this gene is normally found in the membrane of the endoplasmic reticulum (ER). However, upon stress to the ER, the encoded protein is cleaved and the released cytoplasmic transcription factor domain translocates to the nucleus. There it activates the transcription of target genes by binding to box-B elements. [provided by RefSeq, Jun 2013],