

VCC1 Polyclonal Antibody

Cat No: HR1AP12014

For research use only

Overview

Product Name	VCC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	CXCL17 VCC1 UNQ473/PRO842
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>C-X-C motif chemokine ligand 17(CXCL17) Homo sapiens The protein encoded by this gene is a mucosal chemokine that attracts immature dendritic cells and blood monocytes to the lungs. The encoded protein also promotes tumorigenesis through an angiogenic activity. Finally, this protein exhibits strong antimicrobial activity against E. coli, S. aureus, Salmonella, P. aeruginosa, and C. albicans. Two transcript variants, one protein-coding and the other non-protein coding, have been found for this gene. [provided by RefSeq, Dec 2015].</p>