

DSCL1 Polyclonal Antibody

Cat No: HR1AP12013

For research use only

Overview

Product Name	DSCL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	225kDa
GeneID?Human?	DSCAML1 DSCAM2 KIAA1132
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DS cell adhesion molecule like 1(DSCAML1) Homo sapiens The protein encoded by this gene is a member of the Ig superfamily of cell adhesion molecules and is involved in neuronal differentiation. The encoded membrane-bound protein localizes to the cell surface, where it forms aggregates that repel neuronal processes of the same cell type. [provided by RefSeq, Sep 2016],