

# OMD Polyclonal Antibody

Cat No: HR1AP12009

For research use only

## Overview

Product Name	OMD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	OMD SLRR2C UNQ190/PRO216
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be implicated in biomineralization processes. Has a function in binding of osteoblasts via the alpha(V)beta(3)-integrin.,PTM:binds keratan sulfate chains.,similarity:Belongs to the small leucine-rich proteoglycan (SLRP) family. Class II subfamily.,similarity:Contains 12 LRR (leucine-rich) repeats.,subunit:binds the alpha(V)beta(3)-integrin.,tissue specificity:Bone specific.,