

LIN41 Polyclonal Antibody

Cat No: HR1AP12005

For research use only

Overview

Product Name	LIN41 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	TRIM71 LIN41
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tripartite motif containing 71(TRIM71) Homo sapiens The protein encoded by this gene is an E3 ubiquitin-protein ligase that binds with miRNAs and maintains the growth and upkeep of embryonic stem cells. This gene also is involved in the G1-S phase transition of the cell cycle. [provided by RefSeq, Dec 2015],