

VMAC Polyclonal Antibody

Cat No: HR1AP12003

For research use only

Overview

Product Name	VMAC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	VMAC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	subcellular location:Colocalizes with vimentin-type intermediate filaments.,