



SKI Polyclonal Antibody

Cat No: HR1AP11998

For research use only

Overview

Product Name	SKI Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	SKI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SKI proto-oncogene(SKI) Homo sapiens This gene encodes the nuclear protooncogene protein homolog of avian sarcoma viral (v-ski) oncogene. It functions as a repressor of TGF-beta signaling, and may play a role in neural tube development and muscle differentiation. [provided by RefSeq, Oct 2009],