



RGMB Polyclonal Antibody

Cat No: HR1AP11992

For research use only

Overview

Product Name	RGMB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	RGMB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	repulsive guidance molecule family member b(RGMB) Homo sapiens RGMB is a glycosylphosphatidylinositol (GPI)-anchored member of the repulsive guidance molecule family (see RGMA, MIM 607362) and contributes to the patterning of the developing nervous system (Samad et al., 2005 [PubMed 15671031]).[supplied by OMIM, Apr 2009],