

SOX3 Polyclonal Antibody

Cat No: HR1AP11990

For research use only

Overview

Product Name	SOX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	SOX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SRY-box 3(SOX3) Homo sapiens This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. Mutations in this gene have been associated with X-linked mental retardation with growth hormone deficiency. [provided by RefSeq, Jul 2008],</p>