

TMED1 Polyclonal Antibody

Cat No: HR1AP11989

For research use only

Overview

Product Name	TMED1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	TMED1 IL1RL1L IL1RL1LG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transmembrane p24 trafficking protein 1(TMED1) Homo sapiens This gene belongs to the TMED (transmembrane emp24 domain-containing) protein family, which is involved in the vesicular trafficking of proteins. The protein encoded by this gene was identified by its interaction with interleukin 1 receptor-like 1 (IL1RL1) and may play a role in innate immunity. This protein lacks any similarity to other interleukin 1 ligands. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013],