

PSD4 Polyclonal Antibody

Cat No: HR1AP11988

For research use only

Overview

Product Name	PSD4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	116kDa
GeneID?Human?	PSD4 EFA6B TIC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Guanine nucleotide exchange factor for ARF6-like protein. Involved in membrane recycling.,similarity:Contains 1 PH domain.,similarity:Contains 1 SEC7 domain.,subcellular location:Accumulates in dynamic actin-rich membrane ruffles and microvilli-like structures.,tissue specificity:Widely expressed. Highest levels of expression are found in placenta, pancreas, spleen, thymus and peripheral blood.,