

I22R1 Polyclonal Antibody

Cat No: HR1AP11984

For research use only

Overview

Product Name	I22R1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	IL22RA1 IL22R
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>interleukin 22 receptor subunit alpha 1(IL22RA1) Homo sapiens The protein encoded by this gene belongs to the class II cytokine receptor family, and has been shown to be a receptor for interleukin 22 (IL22). IL22 receptor is a protein complex that consists of this protein and interleukin 10 receptor, beta (IL10BR/CRFB4), a subunit also shared by the receptor complex for interleukin 10 (IL10). This gene and interleukin 28 receptor, alpha (IL28RA) form a cytokine receptor gene cluster in the chromosomal region 1p36. [provided by RefSeq, Jul 2008].</p>