

## BECP1 Polyclonal Antibody

Cat No: HR1AP11980

For research use only

### Overview

Product Name	BECP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	BECN1P1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>BECN2 (Beclin 2) is a Protein Coding gene. Among its related pathways are Autophagy Pathway and Kaposi's sarcoma-associated herpesvirus infection. An important paralog of this gene is BECN1. Involved in 2 distinct lysosomal degradation pathways: acts as a regulator of autophagy and as a regulator of G-protein coupled receptors turnover. Regulates degradation in lysosomes of a variety of G-protein coupled receptors via its interaction with GPRASP1/GASP1.</p>