

# ATG9B Polyclonal Antibody

Cat No: HR1AP11979

For research use only

## Overview

Product Name	ATG9B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	ATG9B APG9L2 NOS3AS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>autophagy related 9B(ATG9B) Homo sapiens This gene functions in the regulation of autophagy, a lysosomal degradation pathway. This gene also functions as an antisense transcript in the posttranscriptional regulation of the endothelial nitric oxide synthase 3 gene, which has 3' overlap with this gene on the opposite strand. Mutations in this gene and disruption of the autophagy process have been associated with multiple cancers. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012],</p>