

# ACTBM Polyclonal Antibody

Cat No: HR1AP11966

For research use only

## Overview

Product Name	ACTBM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	POTEKP ACTBL3 FKSG30
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the actin family.,subunit:Interacts with PFN1 and PFDN1. Does not interact with PFN2.,tissue specificity:Expressed in some hepatocellular carcinomas.,