

# SUMO3 Polyclonal Antibody

Cat No: HR1AP11960

For research use only

## Overview

Product Name	SUMO3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	SUMO3 SMT3B SMT3H1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>small ubiquitin-like modifier 3(SUMO3) Homo sapiens This gene encodes a member of the small ubiquitin-related modifier (SUMO) family of eukaryotic proteins. The encoded protein is covalently conjugated to other proteins via a post-translation modification known as sumoylation. Sumoylation may play a role in a wide variety of cellular processes, including nuclear transport, DNA replication and repair, mitosis, transcriptional regulation, and signal transduction. Alternatively spliced transcript variants encoding distinct proteins have been described. [provided by RefSeq, Feb 2014],</p>