

# APOBR Polyclonal Antibody

Cat No: HR1AP11953

For research use only

## Overview

Product Name	APOBR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	119kDa
GeneID?Human?	APOBR APOB48R
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>apolipoprotein B receptor(APOBR) Homo sapiens Apolipoprotein B48 receptor is a macrophage receptor that binds to the apolipoprotein B48 of dietary triglyceride (TG)-rich lipoproteins. This receptor may provide essential lipids, lipid-soluble vitamins and other nutrients to reticuloendothelial cells. If overwhelmed with elevated plasma triglyceride, the apolipoprotein B48 receptor may contribute to foam cell formation, endothelial dysfunction, and atherothrombogenesis. [provided by RefSeq, Jul 2008].</p>