



TSSK1 Polyclonal Antibody

Cat No: HR1AP11946

For research use only

Overview

Product Name	TSSK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	TSSK1B SPOGA1 SPOGA4 STK22A STK22D TSSK1 FKSG81
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	testis specific serine kinase 1B(TSSK1B) Homo sapiens TSSK1 belongs to a family of serine/threonine kinases highly expressed in testis (Hao et al., 2004 [PubMed 15044604]).[supplied by OMIM, Mar 2008],