



TSSK3 Polyclonal Antibody

Cat No: HR1AP11947

For research use only

Overview

Product Name	TSSK3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	TSSK3 SPOGA3 STK22C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	testis specific serine kinase 3(TSSK3) Homo sapiens This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function. [provided by RefSeq, Jul 2008],