

## TRIO Polyclonal Antibody

Cat No: HR1AP11945

For research use only

### Overview

Product Name	TRIO Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	340kDa
GeneID?Human?	TRIO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	trio Rho guanine nucleotide exchange factor(TRIO) Homo sapiens This gene encodes a large protein that functions as a GDP to GTP exchange factor. This protein promotes the reorganization of the actin cytoskeleton, thereby playing a role in cell migration and growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],