

TRIB3 Polyclonal Antibody

Cat No: HR1AP11944

For research use only

Overview

Product Name	TRIB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	TRIB3 C20orf97 NIPK SKIP3 TRB3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tribbles pseudokinase 3(TRIB3) Homo sapiens The protein encoded by this gene is a putative protein kinase that is induced by the transcription factor NF-kappaB. The encoded protein is a negative regulator of NF-kappaB and can also sensitize cells to TNF- and TRAIL-induced apoptosis. In addition, this protein can negatively regulate the cell survival serine-threonine kinase AKT1. Differential promoter usage and alternate splicing result in multiple transcript variants. [provided by RefSeq, Jul 2014],