

TAOK2 Polyclonal Antibody

Cat No: HR1AP11943

For research use only

Overview

Product Name	TAOK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	135kDa
GeneID?Human?	TAOK2 KIAA0881 MAP3K17 PSK PSK1 UNQ2971/PRO7431
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TAO kinase 2(TAOK2) Homo sapiens This gene encodes a serine/threonine protein kinase that is involved in many different processes, including, cell signaling, microtubule organization and stability, and apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011],