

SNRK Polyclonal Antibody

Cat No: HR1AP11941

For research use only

Overview

Product Name	SNRK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	SNRK KIAA0096 SNFRK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SNF related kinase(SNRK) Homo sapiens SNRK is a member of the sucrose nonfermenting (SNF)-related kinase family of serine/threonine kinases (Kertesz et al., 2002 [PubMed 12234663]).[supplied by OMIM, Apr 2009],