

SGK3 Polyclonal Antibody

Cat No: HR1AP11940

For research use only

Overview

Product Name	SGK3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	SGK3 CISK SGKL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serum/glucocorticoid regulated kinase family member 3(SGK3) Homo sapiens This gene is a member of the Ser/Thr protein kinase family and encodes a phosphoprotein with a PX (phox homology) domain. The protein phosphorylates several target proteins and has a role in neutral amino acid transport and activation of potassium and chloride channels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],