

## CDK16 Polyclonal Antibody

Cat No: HR1AP11931

For research use only

### Overview

Product Name	CDK16 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	CDK16 PCTK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cyclin dependent kinase 16(CDK16) Homo sapiens The protein encoded by this gene belongs to the cdc2/cdkx subfamily of the ser/thr family of protein kinases. It may play a role in signal transduction cascades in terminally differentiated cells; in exocytosis; and in transport of secretory cargo from the endoplasmic reticulum. This gene is thought to escape X inactivation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2009].</p>