

NRK Polyclonal Antibody

Cat No: HR1AP11929

For research use only

Overview

Product Name	NRK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	174kDa
GeneID?Human?	NRK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Nik related kinase(NRK) Homo sapiens The mouse ortholog of this gene encodes a protein kinase required for JNK activation. The encoded protein may be involved in the induction of actin polymerization in late embryogenesis.[provided by RefSeq, Jun 2010],