

NEK2 Polyclonal Antibody

Cat No: HR1AP11927

For research use only

Overview

Product Name	NEK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	NEK2 NEK2A NLK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>NIMA related kinase 2(NEK2) Homo sapiens This gene encodes a serine/threonine-protein kinase that is involved in mitotic regulation. This protein is localized to the centrosome, and undetectable during G1 phase, but accumulates progressively throughout the S phase, reaching maximal levels in late G2 phase. Alternatively spliced transcript variants encoding different isoforms with distinct C-termini have been noted for this gene. [provided by RefSeq, Feb 2011],</p>