

## ST17A Polyclonal Antibody

Cat No: HR1AP11923

For research use only

### Overview

Product Name	ST17A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	STK17A DRAK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serine/threonine kinase 17a(STK17A) Homo sapiens This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity. [provided by RefSeq, Jul 2008],