

KC1G3 Polyclonal Antibody

Cat No: HR1AP11920

For research use only

Overview

Product Name	KC1G3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	CSNK1G3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>casein kinase 1 gamma 3(CSNK1G3) Homo sapiens This gene encodes a member of a family of serine/threonine protein kinases that phosphorylate caseins and other acidic proteins. A related protein in the African clawed frog participates in the transmission of Wnt/beta-catenin signaling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],</p>