

## BUB1B Polyclonal Antibody

Cat No: HR1AP11917

For research use only

### Overview

Product Name	BUB1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	115kDa
GeneID?Human?	BUB1B BUBR1 MAD3L SSK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BUB1 mitotic checkpoint serine/threonine kinase B(BUB1B) Homo sapiens This gene encodes a kinase involved in spindle checkpoint function. The protein has been localized to the kinetochore and plays a role in the inhibition of the anaphase-promoting complex/cyclosome (APC/C), delaying the onset of anaphase and ensuring proper chromosome segregation. Impaired spindle checkpoint function has been found in many forms of cancer. [provided by RefSeq, Jul 2008],