

B2L13 Polyclonal Antibody

Cat No: HR1AP11913

For research use only

Overview

Product Name	B2L13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	BCL2L13 MIL1 CD003
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BCL2 like 13(BCL2L13) Homo sapiens This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012].