

NGAP Polyclonal Antibody

Cat No: HR1AP11909

For research use only

Overview

Product Name	NGAP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	RASAL2 NGAP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RAS protein activator like 2(RASAL2) Homo sapiens This gene encodes a protein that contains the GAP-related domain (GRD), a characteristic domain of GTPase-activating proteins (GAPs). GAPs function as activators of Ras superfamily of small GTPases. The protein encoded by this gene is able to complement the defective RasGAP function in a yeast system. Two alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008],</p>