

## RASEF Polyclonal Antibody

Cat No: HR1AP11911

For research use only

### Overview

Product Name	RASEF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	RASEF RAB45
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAS and EF-hand domain containing(RASEF) Homo sapiens This gene is a member of the Rab family of GTPases that are involved in regulation of membrane traffic. The encoded protein contains an N-terminal EF-hand domain, a coiled-coil motif and a C-terminal Rab domain. A potential role as tumor suppressor has been indicated for this gene. [provided by RefSeq, Nov 2012],